

Specimen Information for Zika Virus Testing through Public Health

All requests for testing through Public Health must be approved by calling the Section of Epidemiology (SOE). *Suspected and confirmed Zika infections are reportable to the State — 907-269-8000*

Public health offers and can facilitate testing for Zika Virus. Appropriate collection and labelling of specimens is crucial. Additionally, there are now commercial laboratories which offer testing; providers can coordinate directly with these companies. *Regardless of where testing occurs,* **all** *suspected cases of Zika Virus in Alaska must be reported to the Section of Epidemiology (SOE) within five (5) days at (907) 269-8000.*

Alaska State Public Health Laboratory – PCR and IgM Testing

Time from Symptoms/Exposure	Type of Specimens	Minimum Amounts	Time to Results
<12 Weeks* *Symptomatic individuals should be tested as soon as possible for optimal PCR results	PCR –Serum (required) paired with (optional): Whole Blood (EDTA); Urine; CSF; Amniotic Fluid IgM Serology - Serum	1.5mL of each specimen type submitted for testing	In State: 2-5 business days CDC PRNT testing: 6-8 weeks Note: Positive PCR results are final Positive IgM results are referred to CDC for confirmatory PRNT testing Negative IgM results are final**

**Depending on timing of symptoms/exposure and testing results, follow-up serology testing may be appropriate.

IMPORTANT NOTES

- All specimen types are acceptable for PCR testing up to 12 weeks post symptom onset or post-exposure for asymptomatic pregnant persons. PCR testing provides best results if specimen is collected within 2 weeks of initial symptoms.
- Do not freeze whole blood, keep cold.
- If sender cannot separate blood for serology testing, ask them to spin down serum separator tube prior to sending to the lab, where they can separate and send the serum.
- Please transfer urine, serum, amniotic fluid, or CSF to a plastic tube with screw cap measuring no more than 5 cm tall and approximately 13 mm in diameter (e.g., 1.8 mL cryotube or 2.0 mL microtube).
- All specimens should be kept cold/frozen. Samples may be placed in an insulated container with ice packs.

SPECIMEN SHIPPING

Please ship Monday-Friday during normal business operations. The Alaska State Public Health Laboratory will forward appropriate specimens to CDC laboratory if necessary. Ship specimen and lab requisitions to:

Alaska State Public Health Laboratory 5455 Dr. Martin Luther King Jr. Avenue Anchorage, AK 99507 Ph: 907-334-2100

RESULT REPORTING

ASPHL will notify the ordering provider of positive results by phone; a hard copy of all results will then be sent afterwards.

Some commercial labs are starting to offer testing; SOE will not pay for or ship these specimens. Providers should be aware of caveats associated with these results: <u>http://emergency.cdc.gov/han/han00392.asp</u>

Alaska Section of Epidemiology 3601 C St, Suite 540 Anchorage, AK 99503 907-269-8000 <u>http://dhss.alaska.gov/dph/Epi/Pages/default.aspx</u>